

## Step 1: Download The Shots Installation Files

- Visit *The Group on Immunization Education of the Society of Teachers of Family Medicine* website

<http://www.immunizationed.org/>

- Click on the Shots link



- Click on [Shots 2006 for Pocket PC](#)

### Questions about Shots?

- ◆ [Shots 2006 for Palm OS](#)
- ◆ [Shots 2006 for Pocket PC](#)
- ◆ [Shots 2006 On-Line](#)

- Under Required Microsoft Components – locate and download (save) the file:

[StrongARM Processor for Pocket PC 2003 and later](#)

### Required Microsoft Components

The links below are for the required Microsoft Component do NOT need to install them again. **Pocket PC 2003 sh are running Pocket PC 2003 you will only need SQL !**

- [Microsoft .NET Compact Framework \(MUST BE INSTALLED\)](#) **THE LINKS BELOW ARE FOR SQL SERVER CE ON POCKET PC**
- [StrongARM Processor for Pocket PC 2003 and later](#)

- Under Pocket PC Mobile 5.0 Shots 2006 CAB File – locate and download (save) the file:

[Shots 2006 Mobile 5.0 CAB](#)

### [Shots 2006 MOBILE 5.0 CAB Install](#)

**If you already have SQL Server CE® installed on your PDA, you do not need to install the 5.0 2006 .CAB file below. Otherwise please install the file below under "Required Microsoft Components" BEFORE you come with the Microsoft .NET Compact Framework**


Pocket PC Mobile 5.0 Shots 2006 CAB File

- [Shots 2006 Mobile 5.0 CAB](#)

Need Help?  
  
 CEC  
 679-8870


## Step 2: Copy The Downloaded Files To Your PDA

- Maximize the ActiveSync Window
- Click on the Explore button





- Drag and drop the two files you just downloaded in the Explorer window that opened (this will copy the files onto your PDA)

File names are:  
sqlce.pp3.arm.CAB  
Shots2006.CAB




- If asked to convert the files for the PDA - click on OK



### Step 3: Run the Installations

- Now turn to your PDA
- Click on the Start Menu (see Number 1 on the right)
- Click on Programs (see Number 2 on the right)




- Select File Explorer from your list of Programs



### File Explorer

- Locate the sqlce.pp3.arm.CAB file and tap it to install this file (Note: this file must be installed first)



- Next locate the Shots2006.CAB file and tap it to install this file



Need Help?  
CEC  
679-8870




# Shots 2006 Install (Windows Mobile 5)



- When you begin the installation of each file you will be asked to continue with the installation. To continue - tap yes on your PDA.

**Yes**



- If you have an expansion card inserted in your PDA you will be asked where you'd like to install the files.

Device

will install the files to your device's main memory

\SD Card


will install the files to your expansion card



- You will then see a message letting you know each installation was successful.

Tap ok (located in the top right corner)

**ok**



Need Help?


CEC  
679-8870

# Shots 2006 Install (Windows Mobile 5)

- If you like you may then delete the two files you originally downloaded.


 **sqlce.ppc3.arm**

 **Shots2006**



## Step 4: Begin Using Shots

- Click on the Start menu on your PDA (see Number 1 on the right)
- Click on Programs (see Number 2 on the right)



- Locate Shots 2006 and tap it to begin the program



**Shots 2006**

Need Help?

CEC  
679-8870